



Labrum Repair Anchor


PressLoc™ PK
KnockLoc™ PK
PunchLoc™ PK



KnockLoc™ PK

- Specific material and unique processing allow the distal end of the insert shaft to bend 25 degrees of variable angle
- Flexible entry approaches to achieve the ideal placement position


Product Code	Specification
A302005	Φ 3.0x12mm (2#Blue-White)
A302036	Φ 3.0x12.5mm (2#Blue, bendable)



PressLoc™ PK

- Pre-drilling
- Barbed structure - push in
- Fully threaded design for strong pull-out strength
- Single lead and dual lead options available

Product Code	Specification
A302006	Φ 3.0x13mm (2#Blue-White/White)



PunchLoc™ PK

- Knotless design, eliminating the risk of suture knot impingement
- Easy manipulating, only two steps are required and time saving
- Specially designed for the glenoid cavity to maximize the preservation of glenoid bone
- "Suture-first" technique provides multiple suture options for the surgical procedure
- Visualization and adjustment of suture tension

Product Code	Specification
A302007	Φ 3.0x15.5mm