

AMAZON  **LCP**™


LCP Wrist Fusion Plates 2.7/3.5 mm



System Overview


The AMAZON™ LCP Wrist Fusion Plates 2.7/3.5 mm System consists of Standard Bend Plates, Short Bend Plates and Straight Plates. Standard Bend Plates are for average sized individuals, Short Bend Plates are for small-stature individuals or patients with previous proximal row carpectomy, and Straight Plates are for patients with unusual wrist anatomy or a severely deformed wrist joint. Implants are available in medical grade stainless steel 316L and commercially pure titanium.

Implant Features



Indications

- Wrist arthrodesis and fractures of small bones of carpus
- Posttraumatic arthrosis of the joints of the wrist
- Rheumatoid wrist deformities requiring restoration
- Complex carpal instability
- Postseptic arthritis of the wrist
- Severe unremitting wrist pain related to motion
- Brachial plexus nerve palsies
- Tumor resection
- Spastic deformities



Implant Information



LCP Wrist Fusion Plate 2.7/3.5 mm, Standard Bend

REF (S.S)	REF (TIT)	HOLES	LENGTH
MOI 15716118	MOI 35716118	8	118 mm

For average sized individuals



LCP Wrist Fusion Plate 2.7/3.5 mm, Short Bend

REF (S.S)	REF (TIT)	HOLES	LENGTH
MOI 15717118	MOI 35717118	8	118 mm

For small-stature individuals or patients with previous proximal row carpectomy



LCP Wrist Fusion Plate 2.7/3.5 mm, Straight


REF (S.S)	REF (TIT)	HOLES	LENGTH
MOI 15718118	MOI 35718112	9	112 mm

For patients with unusual wrist anatomy or a severely deformed wrist joint



2.7 mm Cortex Screws

REF (S.S)	REF (TIT)	LENGTH
MOI 12056010	MOI 32056010	10 mm
MOI 12056012	MOI 32056012	12 mm
MOI 12056014	MOI 32056014	14 mm
MOI 12056016	MOI 32056016	16 mm
MOI 12056018	MOI 32056018	18 mm
MOI 12056020	MOI 32056020	20 mm
MOI 12056022	MOI 32056022	22 mm
MOI 12056024	MOI 32056024	24 mm
MOI 12056026	MOI 32056026	26 mm



2.7 mm Locking Screws

REF (S.S)	REF (TIT)	LENGTH
MOI 15224010	MOI 35224010	10 mm
MOI 15224012	MOI 35224012	12 mm
MOI 15224014	MOI 35224014	14 mm
MOI 15224016	MOI 35224016	16 mm
MOI 15224018	MOI 35224018	18 mm
MOI 15224020	MOI 35224020	20 mm
MOI 15224022	MOI 35224022	22 mm
MOI 15224024	MOI 35224024	24 mm
MOI 15224026	MOI 35224026	26 mm



3.5 mm Locking Screws

REF (S.S)	REF (TIT)	LENGTH
MOI 15228012	MOI 35228012	12 mm
MOI 15228014	MOI 35228014	14 mm
MOI 15228016	MOI 35228016	16 mm
MOI 15228018	MOI 35228018	18 mm
MOI 15228020	MOI 35228020	20 mm
MOI 15228022	MOI 35228022	22 mm
MOI 15228024	MOI 35228024	24 mm
MOI 15228026	MOI 35228026	26 mm
MOI 15228028	MOI 35228028	28 mm
MOI 15228030	MOI 35228030	30 mm



3.5 mm Cortex Screws

REF (S.S)	REF (TIT)	LENGTH
MOI 13500012	MOI 33500012	12 mm
MOI 13500014	MOI 33500014	14 mm
MOI 13500016	MOI 33500016	16 mm
MOI 13500018	MOI 33500018	18 mm
MOI 13500020	MOI 33500020	20 mm
MOI 13500022	MOI 33500022	22 mm
MOI 13500024	MOI 33500024	24 mm
MOI 13500026	MOI 33500026	26 mm
MOI 13500028	MOI 33500028	28 mm
MOI 13500030	MOI 33500030	30 mm



| TRAUMA

USA OFFICE

2711 S Ocean Dr.
Unit 2704
Hollywood, FL 33019

PUERTO RICO OFFICE

1660 Santa Ana Street
Esq. Fidalgo Diaz,
San Juan, Puerto Rico 00909

E-mail: info@madisonorthoinc.com
Website: www.madisonorthoinc.com

