

CANNA

Cannulated Screws System



The CANNA™ Cannulated Screw System is intended for fixation of fractures, fusions, and osteotomies of large and small bones appropriate for the size of the implants. All screws are cannulated in order to be used over a guide wire, and each screw is either partially or fully threaded using a cancellous thread form. Partially threaded screws may be used to lag one bone fragment to another, where the bone fragment is captured by the threads of the screw and pulled toward the near cortex fragment on the head side of the screw. Fully threaded screws are intended to be used to stabilize fractures with little to no compression across the fracture.

Indications for Use

■ 2.5 mm Cannulated Screws

- Intra-articular fractures of the carpals, metacarpals, tarsals and metatarsals
- Fixation of small bone fragments
- Bunionectomies and osteotomies
- Arthrodeses of small joints



■ 3.0 mm Cannulated Screws

- Fixation of fractures of the forearm, hand and foot, e.g.:
- Fractures and arthrodeses of the carpals and metacarpals
 - Fractures of the distal radius and radial head
 - Metatarsal fractures



■ 3.5 mm Cannulated Screws

- Fixation of fractures with small fragments, e.g.:
- Wrist fractures
 - Metacarpal and metatarsal fractures and fixation in metacarpal and metatarsal osteotomies
 - Tarsal fractures
 - Transcondylar humeral fractures in children




■ 4.0 mm Cannulated Screws

- Fixation of fractures with medium fragments, e.g.:
- Tarsal and metatarsal fractures and fixation in metatarsal and phalangeal osteotomies
 - Tarsometatarsal and metatarsophalangeal arthrodeses
 - Ligament fixations
 - Hallux valgus corrections



■ 4.5 mm/5.5 mm Cannulated Screws

- Fixation of fractures with medium fragments, e.g.:
- Malleolar fractures
 - Pilon tibial fractures
 - Fractures of the calcaneus and talus
 - Tibial plateau fractures
 - Carpal and tarsal arthrodeses



■ 6.5 mm/7.5 mm Cannulated Screws

- Fixation of fractures with large fragments, e.g.:
- Femoral neck fractures
 - Intercondylar femoral fractures
 - Epiphyseolysis of the femoral head
 - Ankle arthrodeses
 - Iliosacral dislocations

IMPLANT INFORMATION

Washer



REF (TIT)	OUTER DIA	TYPE
MOI 36399065	6.5	For 2.5 mm Screws
MOI 36399009	9.0	For 3.0 / 3.5 / 4.0 mm Screws
MOI 36399010	10.0	For 4.5 / 5.5 mm Screws
MOI 36399013	13.0	For 6.5 / 7.5 mm Screws

2.5 mm Cannulated Screw, Short Thread



REF (TIT)	SCREW LENGTH	THREAD LENGTH
MOI 36402010	10 mm	3 mm
MOI 36402012	12 mm	3 mm
MOI 36402014	14 mm	3 mm
MOI 36402016	16 mm	4 mm
MOI 36402018	18 mm	4 mm
MOI 36402020	20 mm	6 mm
MOI 36402022	22 mm	6 mm
MOI 36402024	24 mm	6 mm
MOI 36402026	26 mm	6 mm
MOI 36402028	28 mm	6 mm
MOI 36402030	30 mm	6 mm

2.5 mm Cannulated Screw, Long Thread



REF (TIT)	SCREW LENGTH	THREAD LENGTH
MOI 36401010	10 mm	5 mm
MOI 36401012	12 mm	5 mm
MOI 36401014	14 mm	5 mm
MOI 36401016	16 mm	5 mm
MOI 36401018	18 mm	5 mm
MOI 36401020	20 mm	12 mm
MOI 36401022	22 mm	12 mm
MOI 36401024	24 mm	12 mm
MOI 36401026	26 mm	12 mm
MOI 36401028	28 mm	12 mm
MOI 36401030	30 mm	12 mm

3.0 mm Cannulated Screw, Short Thread



REF (TIT)	SCREW LENGTH	THREAD LENGTH
MOI 36404010	10 mm	4 mm
MOI 36404012	12 mm	4 mm
MOI 36404014	14 mm	4 mm
MOI 36404016	16 mm	4 mm
MOI 36404018	18 mm	6 mm
MOI 36404020	20 mm	6 mm
MOI 36404022	22 mm	6 mm
MOI 36404024	24 mm	6 mm
MOI 36404026	26 mm	6 mm
MOI 36404028	28 mm	8 mm
MOI 36404030	30 mm	8 mm
MOI 36404032	32 mm	8 mm
MOI 36404034	34 mm	8 mm
MOI 36404036	36 mm	8 mm
MOI 36404038	38 mm	8 mm
MOI 36404040	40 mm	8 mm

3.0 mm Cannulated Screw, Long Thread



REF (TIT)	SCREW LENGTH	THREAD LENGTH
MOI 36403014	14 mm	8 mm
MOI 36403016	16 mm	8 mm
MOI 36403018	18 mm	8 mm
MOI 36403020	20 mm	8 mm
MOI 36403022	22 mm	8 mm
MOI 36403024	24 mm	8 mm
MOI 36403026	26 mm	16 mm
MOI 36403028	28 mm	16 mm
MOI 36403030	30 mm	16 mm
MOI 36403032	32 mm	16 mm
MOI 36403034	34 mm	16 mm
MOI 36403036	36 mm	16 mm
MOI 36403038	38 mm	16 mm
MOI 36403040	40 mm	16 mm

IMPLANT INFORMATION

3.5 mm Cannulated Screw, Short Thread



REF (TIT)	SCREW LENGTH	THREAD LENGTH
MOI 36406012	12 mm	4 mm
MOI 36406014	14 mm	4 mm
MOI 36406016	16 mm	4 mm
MOI 36406018	18 mm	4 mm
MOI 36406020	20 mm	6 mm
MOI 36406022	22 mm	6 mm
MOI 36406024	24 mm	6 mm
MOI 36406026	26 mm	6 mm
MOI 36406028	28 mm	6 mm
MOI 36406030	30 mm	8 mm
MOI 36406032	32 mm	8 mm
MOI 36406034	34 mm	8 mm
MOI 36406036	36 mm	8 mm
MOI 36406038	38 mm	8 mm
MOI 36406040	40 mm	8 mm

3.5 mm Cannulated Screw, Full Thread



REF (TIT)	SCREW LENGTH	THREAD LENGTH
MOI 36405012	12 mm	7 mm
MOI 36405014	14 mm	9 mm
MOI 36405016	16 mm	11 mm
MOI 36405018	18 mm	13 mm
MOI 36405020	20 mm	15 mm
MOI 36405022	22 mm	17 mm
MOI 36405024	24 mm	19 mm
MOI 36405026	26 mm	21 mm
MOI 36405028	28 mm	23 mm
MOI 36405030	30 mm	25 mm
MOI 36405032	32 mm	27 mm
MOI 36405034	34 mm	29 mm
MOI 36405035	35 mm	30 mm
MOI 36405036	36 mm	31 mm
MOI 36405038	38 mm	33 mm
MOI 36405040	40 mm	35 mm
MOI 36405045	45 mm	40 mm
MOI 36405050	50 mm	45 mm

3.5 mm Cannulated Screw, Long Thread



REF (TIT)	SCREW LENGTH	THREAD LENGTH
MOI 36418014	14 mm	8 mm
MOI 36418016	16 mm	8 mm
MOI 36418018	18 mm	8 mm
MOI 36418020	20 mm	8 mm
MOI 36418022	22 mm	8 mm
MOI 36418024	24 mm	8 mm
MOI 36418026	26 mm	16 mm
MOI 36418028	28 mm	16 mm
MOI 36418030	30 mm	16 mm
MOI 36418032	32 mm	16 mm
MOI 36418034	34 mm	16 mm
MOI 36418035	35 mm	16 mm
MOI 36418036	36 mm	16 mm
MOI 36418038	38 mm	16 mm
MOI 36418040	40 mm	16 mm
MOI 36418045	45 mm	16 mm
MOI 36418050	50 mm	16 mm

4.0 mm Cannulated Screw, Short Thread



REF (TIT)	SCREW LENGTH	THREAD LENGTH
MOI 36408018	18 mm	6 mm
MOI 36408020	20 mm	6 mm
MOI 36408022	22 mm	6 mm
MOI 36408024	24 mm	6 mm
MOI 36408026	26 mm	6 mm
MOI 36408028	28 mm	6 mm
MOI 36408030	30 mm	8 mm
MOI 36408032	32 mm	8 mm
MOI 36408034	34 mm	8 mm
MOI 36408035	35 mm	8 mm
MOI 36408036	36 mm	8 mm
MOI 36408038	38 mm	8 mm
MOI 36408040	40 mm	10 mm
MOI 36408042	42 mm	10 mm
MOI 36408044	44 mm	10 mm
MOI 36408046	46 mm	10 mm
MOI 36408048	48 mm	10 mm
MOI 36408050	50 mm	10 mm

IMPLANT INFORMATION

4.0 mm Cannulated Screw, Long Thread



REF (TIT)	SCREW LENGTH	THREAD LENGTH
MOI 36407018	18 mm	10 mm
MOI 36407020	20 mm	10 mm
MOI 36407022	22 mm	10 mm
MOI 36407024	24 mm	10 mm
MOI 36407026	26 mm	10 mm
MOI 36407028	28 mm	10 mm
MOI 36407030	30 mm	10 mm
MOI 36407032	32 mm	16 mm
MOI 36407034	34 mm	16 mm
MOI 36407035	35 mm	16 mm
MOI 36407036	36 mm	16 mm
MOI 36407038	38 mm	16 mm
MOI 36407040	40 mm	16 mm
MOI 36407042	42 mm	16 mm
MOI 36407044	44 mm	16 mm
MOI 36407046	46 mm	16 mm
MOI 36407048	48 mm	16 mm
MOI 36407050	50 mm	16 mm

4.5 mm Cannulated Screw, Short Thread



REF (TIT)	SCREW LENGTH	THREAD LENGTH
MOI 36409020	20 mm	7 mm
MOI 36409022	22 mm	7 mm
MOI 36409024	24 mm	8 mm
MOI 36409026	26 mm	9 mm
MOI 36409028	28 mm	9 mm
MOI 36409030	30 mm	10 mm
MOI 36409032	32 mm	11 mm
MOI 36409034	34 mm	11 mm
MOI 36409036	36 mm	12 mm
MOI 36409038	38 mm	13 mm
MOI 36409040	40 mm	13 mm
MOI 36409042	42 mm	14 mm
MOI 36409044	44 mm	15 mm
MOI 36409046	46 mm	15 mm
MOI 36409048	48 mm	16 mm
MOI 36409050	50 mm	16 mm
MOI 36409052	52 mm	16 mm
MOI 36409054	54 mm	16 mm
MOI 36409056	56 mm	16 mm
MOI 36409058	58 mm	16 mm
MOI 36409060	60 mm	16 mm

4.5 mm Cannulated Screw, Long Thread



REF (TIT)	SCREW LENGTH	THREAD LENGTH
MOI 36419020	20 mm	12 mm
MOI 36419022	22 mm	12 mm
MOI 36419024	24 mm	12 mm
MOI 36419026	26 mm	12 mm
MOI 36419028	28 mm	12 mm
MOI 36419030	30 mm	12 mm
MOI 36419032	32 mm	12 mm
MOI 36419034	34 mm	12 mm
MOI 36419036	36 mm	12 mm
MOI 36419038	38 mm	12 mm
MOI 36419040	40 mm	12 mm
MOI 36419042	42 mm	16 mm
MOI 36419044	44 mm	16 mm
MOI 36419046	46 mm	16 mm
MOI 36419048	48 mm	16 mm
MOI 36419050	50 mm	16 mm
MOI 36419052	52 mm	16 mm
MOI 36419054	54 mm	16 mm
MOI 36419056	56 mm	16 mm
MOI 36419058	58 mm	16 mm
MOI 36419060	60 mm	16 mm

4.5 mm Cannulated Screw, Full Thread



REF (TIT)	SCREW LENGTH	THREAD LENGTH
MOI 36410020	20 mm	15 mm
MOI 36410022	22 mm	17 mm
MOI 36410024	24 mm	19 mm
MOI 36410026	26 mm	21 mm
MOI 36410028	28 mm	23 mm
MOI 36410030	30 mm	25 mm
MOI 36410032	32 mm	27 mm
MOI 36410034	34 mm	29 mm
MOI 36410036	36 mm	31 mm
MOI 36410038	38 mm	33 mm
MOI 36410040	40 mm	35 mm
MOI 36410042	42 mm	37 mm
MOI 36410044	44 mm	39 mm
MOI 36410046	46 mm	41 mm
MOI 36410048	48 mm	43 mm
MOI 36410050	50 mm	45 mm
MOI 36410052	52 mm	47 mm
MOI 36410054	54 mm	49 mm
MOI 36410056	56 mm	51 mm
MOI 36410058	58 mm	53 mm
MOI 36410060	60 mm	55 mm

IMPLANT INFORMATION

5.5 mm Cannulated Screw, Short Thread



REF (TIT)	SCREW LENGTH	THREAD LENGTH
MOI 36412020	20 mm	7 mm
MOI 36412022	22 mm	7 mm
MOI 36412024	24 mm	8 mm
MOI 36412026	26 mm	9 mm
MOI 36412028	28 mm	9 mm
MOI 36412030	30 mm	10 mm
MOI 36412032	32 mm	11 mm
MOI 36412034	34 mm	11 mm
MOI 36412036	36 mm	12 mm
MOI 36412038	38 mm	13 mm
MOI 36412040	40 mm	13 mm
MOI 36412042	42 mm	14 mm
MOI 36412044	44 mm	15 mm
MOI 36412046	46 mm	15 mm
MOI 36412048	48 mm	16 mm
MOI 36412050	50 mm	16 mm
MOI 36412052	52 mm	16 mm
MOI 36412054	54 mm	16 mm
MOI 36412056	56 mm	16 mm
MOI 36412058	58 mm	16 mm
MOI 36412060	60 mm	16 mm
MOI 36412065	65 mm	16 mm
MOI 36412070	70 mm	16 mm

5.5 mm Cannulated Screw, Long Thread



REF (TIT)	SCREW LENGTH	THREAD LENGTH
MOI 36411020	20 mm	12 mm
MOI 36411022	22 mm	12 mm
MOI 36411024	24 mm	12 mm
MOI 36411026	26 mm	12 mm
MOI 36411028	28 mm	14 mm
MOI 36411030	30 mm	14 mm
MOI 36411032	32 mm	14 mm
MOI 36411034	34 mm	14 mm
MOI 36411036	36 mm	14 mm
MOI 36411038	38 mm	14 mm
MOI 36411040	40 mm	32 mm
MOI 36411042	42 mm	32 mm
MOI 36411044	44 mm	32 mm
MOI 36411046	46 mm	32 mm
MOI 36411048	48 mm	32 mm
MOI 36411050	50 mm	32 mm
MOI 36411052	52 mm	32 mm
MOI 36411054	54 mm	32 mm
MOI 36411056	56 mm	32 mm
MOI 36411058	58 mm	32 mm
MOI 36411060	60 mm	32 mm
MOI 36411065	65 mm	32 mm
MOI 36411070	70 mm	32 mm

6.5 mm Cannulated Screw, Short Thread



REF (TIT)	SCREW LENGTH	THREAD LENGTH
MOI 36414030	30 mm	16 mm
MOI 36414035	35 mm	16 mm
MOI 36414040	40 mm	16 mm
MOI 36414045	45 mm	16 mm
MOI 36414050	50 mm	16 mm
MOI 36414055	55 mm	16 mm
MOI 36414060	60 mm	16 mm
MOI 36414065	65 mm	16 mm
MOI 36414070	70 mm	16 mm
MOI 36414075	75 mm	16 mm
MOI 36414080	80 mm	16 mm
MOI 36414085	85 mm	16 mm
MOI 36414090	90 mm	16 mm
MOI 36414095	95 mm	16 mm
MOI 36414100	100 mm	16 mm
MOI 36414105	105 mm	16 mm
MOI 36414110	110 mm	16 mm
MOI 36414115	115 mm	16 mm
MOI 36414120	120 mm	16 mm
MOI 36414125	125 mm	16 mm

6.5 mm Cannulated Screw, Long Thread



REF (TIT)	SCREW LENGTH	THREAD LENGTH
MOI 36413045	45 mm	32 mm
MOI 36413050	50 mm	32 mm
MOI 36413055	55 mm	32 mm
MOI 36413060	60 mm	32 mm
MOI 36413065	65 mm	32 mm
MOI 36413070	70 mm	32 mm
MOI 36413075	75 mm	32 mm
MOI 36413080	80 mm	32 mm
MOI 36413085	85 mm	32 mm
MOI 36413090	90 mm	32 mm
MOI 36413095	95 mm	32 mm
MOI 36413100	100 mm	32 mm
MOI 36413105	105 mm	32 mm
MOI 36413110	110 mm	32 mm
MOI 36413115	115 mm	32 mm
MOI 36413120	120 mm	32 mm
MOI 36413125	125 mm	32 mm

IMPLANT INFORMATION

7.0 mm Cannulated Screw, Short Thread



REF (TIT)	SCREW LENGTH	THREAD LENGTH
MOI 35750030	30 mm	16 mm
MOI 35750035	35 mm	16 mm
MOI 35750040	40 mm	16 mm
MOI 35750045	45 mm	16 mm
MOI 35750050	50 mm	16 mm
MOI 35750055	55 mm	16 mm
MOI 35750060	60 mm	16 mm
MOI 35750065	65 mm	16 mm
MOI 35750070	70 mm	16 mm
MOI 35750075	75 mm	16 mm
MOI 35750080	80 mm	16 mm
MOI 35750085	85 mm	16 mm
MOI 35750090	90 mm	16 mm
MOI 35750095	95 mm	16 mm
MOI 35750100	100 mm	16 mm
MOI 35750105	105 mm	16 mm
MOI 35750110	110 mm	16 mm
MOI 35750115	115 mm	16 mm
MOI 35750120	120 mm	16 mm
MOI 35750125	125 mm	16 mm
MOI 35750130	130 mm	16 mm
MOI 35750135	135 mm	16 mm

7.0 mm Cannulated Screw, Long Thread



REF (TIT)	SCREW LENGTH	THREAD LENGTH
MOI 35751045	45 mm	32 mm
MOI 35751050	50 mm	32 mm
MOI 35751055	55 mm	32 mm
MOI 35751060	60 mm	32 mm
MOI 35751065	65 mm	32 mm
MOI 35751070	70 mm	32 mm
MOI 35751075	75 mm	32 mm
MOI 35751080	80 mm	32 mm
MOI 35751085	85 mm	32 mm
MOI 35751090	90 mm	32 mm
MOI 35751095	95 mm	32 mm
MOI 35751100	100 mm	32 mm
MOI 35751105	105 mm	32 mm
MOI 35751110	110 mm	32 mm
MOI 35751115	115 mm	32 mm
MOI 35751120	120 mm	32 mm
MOI 35751125	125 mm	32 mm
MOI 35751130	130 mm	32 mm
MOI 35751135	135 mm	32 mm

7.5 mm Cannulated Screw, Short Thread



REF (TIT)	SCREW LENGTH	THREAD LENGTH
MOI 36416030	30 mm	16 mm
MOI 36416035	35 mm	16 mm
MOI 36416040	40 mm	16 mm
MOI 36416045	45 mm	16 mm
MOI 36416050	50 mm	16 mm
MOI 36416055	55 mm	16 mm
MOI 36416060	60 mm	16 mm
MOI 36416065	65 mm	16 mm
MOI 36416070	70 mm	16 mm
MOI 36416075	75 mm	16 mm
MOI 36416080	80 mm	16 mm
MOI 36416085	85 mm	16 mm
MOI 36416090	90 mm	16 mm
MOI 36416095	95 mm	16 mm
MOI 36416100	100 mm	16 mm
MOI 36416105	105 mm	16 mm
MOI 36416110	110 mm	16 mm
MOI 36416115	115 mm	16 mm
MOI 36416120	120 mm	16 mm
MOI 36416125	125 mm	16 mm
MOI 36416130	130 mm	16 mm
MOI 36416135	135 mm	16 mm

7.5 mm Cannulated Screw, Long Thread



REF (TIT)	SCREW LENGTH	THREAD LENGTH
MOI 36415045	45 mm	32 mm
MOI 36415050	50 mm	32 mm
MOI 36415055	55 mm	32 mm
MOI 36415060	60 mm	32 mm
MOI 36415065	65 mm	32 mm
MOI 36415070	70 mm	32 mm
MOI 36415075	75 mm	32 mm
MOI 36415080	80 mm	32 mm
MOI 36415085	85 mm	32 mm
MOI 36415090	90 mm	32 mm
MOI 36415095	95 mm	32 mm
MOI 36415100	100 mm	32 mm
MOI 36415105	105 mm	32 mm
MOI 36415110	110 mm	32 mm
MOI 36415115	115 mm	32 mm
MOI 36415120	120 mm	32 mm
MOI 36415125	125 mm	32 mm
MOI 36415130	130 mm	32 mm
MOI 36415135	135 mm	32 mm



| TRAUMA

USA OFFICE

2711 S Ocean Dr.
Unit 2704
Hollywood, FL 33019

PUERTO RICO OFFICE

1660 Santa Ana Street
Esq. Fidalgo Diaz,
San Juan, Puerto Rico 00909

E-mail: info@madisonorthoinc.com
Website: www.madisonorthoinc.com

