


ROCKY™

Tibia Nailing System





INTRODUCTION

ROCKY™ Tibia Nailing System is suitable for all fractures extending from 7-8 cm distal to the tibial plateau, to within 5.5 cm of the distal articular surface, provided that the epiphyses are closed. ROCKY™ is a set of intramedullary nails which offers several advantages over existing systems. A major advantage of the system is the ability to insert both proximal and distal locking screws accurately and quickly without the use of X-rays, using an external mechanical targeting device.

Features and Benefits

- Minimal invasive technique and biomechanical stability allow early functional exercise.
- Accurate and fast mechanical aiming device for proximal and distal locking.
- A locking screw of a smooth shank configuration is much stronger for a given diameter than a fully threaded screw.
- All of the implants are available in Titanium Alloy and Stainless Steel.
- The nails and locking screws are color coded for easy identification.



Indications

The ROCKY™ Tibia Nailing System is indicated for: open or closed tibia shaft fractures with a very proximal and very distal extent in which locking screw fixation can be obtained, multi-fragment fractures, segmental fractures, proximal or distal non-unions, proximal or distal mal-unions, pseudarthrosis, corrective osteotomies, pathologic and impending pathologic fractures, tumor resections and comminuted fractures with or without bone loss.



IMPLANT INFORMATION



ROCKY™ Tibia Nail, Cannulated

REF (SS)	REF (TIT)	DIA (mm)	LEN (mm)
MOI 11301240	MOI 31301240	8.0	240
MOI 11301260	MOI 31301260	8.0	260
MOI 11301280	MOI 31301280	8.0	280
MOI 11301300	MOI 31301300	8.0	300
MOI 11301320	MOI 31301320	8.0	320
MOI 11301340	MOI 31301340	8.0	340
MOI 11301360	MOI 31301360	8.0	360
MOI 11302240	MOI 31302240	9.0	240
MOI 11302260	MOI 31302260	9.0	260
MOI 11302280	MOI 31302280	9.0	280
MOI 11302300	MOI 31302300	9.0	300
MOI 11302320	MOI 31302320	9.0	320
MOI 11302340	MOI 31302340	9.0	340
MOI 11302360	MOI 31302360	9.0	360
MOI 11303240	MOI 31303240	10.0	240
MOI 11303260	MOI 31303260	10.0	260
MOI 11303280	MOI 31303280	10.0	280
MOI 11303300	MOI 31303300	10.0	300
MOI 11303320	MOI 31303320	10.0	320
MOI 11303340	MOI 31303340	10.0	340
MOI 11303360	MOI 31303360	10.0	360

ROCKY™ Tibia Nail, Solid

REF (SS)	REF (TIT)	DIA (mm)	LEN (mm)
MOI 11001240	MOI 31001240	8.0	240
MOI 11001260	MOI 31001260	8.0	260
MOI 11001280	MOI 31001280	8.0	280
MOI 11001300	MOI 31001300	8.0	300
MOI 11001320	MOI 31001320	8.0	320
MOI 11001340	MOI 31001340	8.0	340
MOI 11001360	MOI 31001360	8.0	360
MOI 11002240	MOI 31002240	9.0	240
MOI 11002260	MOI 31002260	9.0	260
MOI 11002280	MOI 31002280	9.0	280
MOI 11002300	MOI 31002300	9.0	300
MOI 11002320	MOI 31002320	9.0	320
MOI 11002340	MOI 31002340	9.0	340
MOI 11002360	MOI 31002360	9.0	360
MOI 11003240	MOI 31003240	10.0	240
MOI 11003260	MOI 31003260	10.0	260
MOI 11003280	MOI 31003280	10.0	280
MOI 11003300	MOI 31003300	10.0	300
MOI 11003320	MOI 31003320	10.0	320
MOI 11003340	MOI 31003340	10.0	340
MOI 11003360	MOI 31003360	10.0	360



Locking Screw 4.0 mm

REF (SS)	REF (TIT)	DIA (mm)	LEN (mm)
MOI 11501030	MOI 31501030	4.0	30
MOI 11501035	MOI 31501035	4.0	35
MOI 11501040	MOI 31501040	4.0	40
MOI 11501045	MOI 31501045	4.0	45
MOI 11501050	MOI 31501050	4.0	50
MOI 11501055	MOI 31501055	4.0	55
MOI 11501060	MOI 31501060	4.0	60
MOI 11501065	MOI 31501065	4.0	65
MOI 11501070	MOI 31501070	4.0	70



End Cap, for ROCKY™ Tibia Nail

REF (SS)	REF (TIT)	SIZE (mm)
MOI 11701000	MOI 31701000	2



USA OFFICE

2404 Hollywood Blvd.
Hollywood, FL 33020

PUERTO RICO OFFICE

1660 Santa Ana Street
Esq. Fidalgo Diaz,
San Juan, Puerto Rico 00909

E-mail: info@madisonorthoinc.com
Website: www.madisonorthoinc.com

