





Cranio-maxillofacial Repair and Reconstruction System



Catalogue | *CMF*



Cranio-maxillofacial Repair and Reconstruction System



2.4
LOCK

Thickness 2.0 - 2.5 mm


Mandible



LOCK Mandible Plate 2.4, straight

Titanium	Holes	Curve Length
MOI 22048006	6	50 mm
MOI 22048012	12	98 mm
MOI 22048020	20	160 mm


Thickness: 2.0 mm Width: 8.0 mm



LOCK Mandible Plate 2.4, angled

Titanium	Holes	L / R
MOI 22049621	6 + 21	Left
MOI 22050621	6 + 21	Right


Thickness: 2.5 mm Width: 8.0mm Curve Length: 214 mm



LOCK Mandible Plate 2.4, double angled

Titanium	Holes	Curve Length
MOI 22051420	4 + 20 + 4	111 mm
MOI 22051522	5 + 22 + 5	119 mm
MOI 22051624	6 + 24 + 6	133 mm

Thickness: 2.5 mm Width: 8.0 mm



LOCK Mandible Plate 2.4, crescent - shaped

Titanium	Holes	Curve Length
MOI 22052022	2 + 2	50 mm
MOI 22052033	3 + 3	70 mm

Thickness: 2.0 mm Width: 8.0 mm



LOCK Mandible Plate 2.4, arc - shaped


Titanium	Holes	Curve Length
MOI 22053033	3 + 3	58 mm
MOI 22053044	4 + 4	79 mm

Thickness: 2.0 mm Width: 8.0 mm




LOCK Screw 2.4 mm, cruciform, self - tapping

TAN	Length (mm)
MOI 32054005	5
MOI 32054006	6
MOI 32054007	7
MOI 32054008	8
MOI 32054010	10
MOI 32054012	12
MOI 32054014	14
MOI 32054016	16
MOI 32054018	18






Cranio-maxillofacial Repair and Reconstruction System




CMF Instrument Sets




Compact CMF Instrument Set

	REF CODE	DESCRIPTION	UNITS
1	MOI102000	Depth Gauge	1
2	MOI1102001	Screwdriver Shaft, cruciform	1
3	MOI1102002	Screw Extractor Ø1.5	1
4	MOI1102003	Screw Extractor Ø2.0	1
5	MOI1102004	Drill Bit Ø1.1 mm	2
6	MOI1102005	Drill Guide Ø1.1 / Ø1.6	1
7	MOI1102006	Drill Bit Ø1.6 mm	2
8	MOI1102007	Handle, with Mini Quick Coupling	1

Compact CMF Instrument Set, for Reconstruction



REF CODE	DESCRIPTION	UNITS
1	MOI1101000 Plate Cutter	2
2	MOI1101001 Bending Iron	2
3	MOI1101002 Handle, with Mini Quick Coupling	1
4	MOI1101003 Depth Gauge	1
5	MOI1101004 Drill Sleeve, with thread	3
6	MOI1101005 Screwdriver Shaft, cruciform	1
7	MOI1101006 Drill Bit Ø1.7 mm	2
8	MOI1101007 Drill Bit Ø1.9 mm	2
9	MOI1101008 Screw Container	1



REF CODE	DESCRIPTION	UNITS
10	MOI1101009 Drill Guide Ø1.6 / Ø1.9	1
11	MOI1101010 Combination Bending Pliers	1
12	MOI1101011 Drill Bit Ø2.0 mm	1
13	MOI1101012 Tap	2



TRAUMA

USA OFFICE

4000 Island Blvd, Suite 2704
Aventura, Florida 33160

PUERTO RICO OFFICE

1660 Santa Ana Street
Esq. Fidalgo Diaz,
San Juan, Puerto Rico 00909

E-mail: info@madisonorthoinc.com
Website: www.madisonorthoinc.com

