


URANUS™


Hammertoe Fixation System



INTRODUCTION

URANUS™ Hammertoe Fixation System is a simple and effective style of hammertoe implant technology. This technique describes the “V”-shaped blade-screw implant that seeks to provide optimal fixation in the proximal and middle phalanx using an “open” technique. Available in 5 sizes, the URANUS™ system provides sizing versatility for a variety of indications and anatomical conditions.

The stainless steel construction of the implants balance strength and workability through the unique design of the one-piece implant. This hybrid design provides secure, multi-axial fixation of hammertoe correction.



Indications for Use

The URANUS™ Hammertoe Implant System is indicated for the fixation of osteotomies and reconstruction of the lesser toes following correction procedures for hammertoe, claw toe and mallet toe.




IMPLANT INFORMATION

URANUS™ Hammertoe Fixation Screws

REF (S.S)	SIZES	ANGLE	
MOI 16021009	2.8 X 9 MM	0°	
MOI 16022009	2.8 X 9 MM	10°	
MOI 16023013	2.0 X 13 MM	0°	
MOI 16024013	2.0 X 13 MM	10°	
MOI 16025016	2.4 X 16 MM	0°	
MOI 16026016	2.4 X 16 MM	10°	
MOI 16027016	3.2 X 16 MM	0°	
MOI 16028016	3.2 X 16 MM	10°	
MOI 16029016	4.0 X 16 MM	0°	
MOI 16030016	4.0 X 16 MM	10°	

INSTRUMENT INFORMATION



REF	DESCRIPTION	UNITS
MOI 1700001U	AO Driver Handle	1
MOI 1700002U	Driver Shaft 0°	1
MOI 1700003U	Driver Shaft 10°	1
MOI 1700004U	K-Wire 1.2 mm	1
MOI 1700005U	Broach-X-Large (3.0 mm and 4.0 mm)	1
MOI 1700006U	Broach-Standard (2.0mm m and 2.4 mm)	1
MOI 1700007U	Broach-X-Small (2.8 mm)	1
MOI 1700008U	6.5 Planar	1
MOI 1700009U	7.5 Planar	1
MOI 1700000U	Tray and Container	1



USA OFFICE

2404 Hollywood Blvd.
Hollywood, FL 33020

PUERTO RICO OFFICE

1660 Santa Ana Street
Esq. Fidalgo Diaz,
San Juan, Puerto Rico 00909

E-mail: info@madisonorthoinc.com
Website: www.madisonorthoinc.com

