

ACCK Revision Knee System

- ACCK is intended for patients who, in the surgeon's judgment, require additional prosthetic stabilization due to inadequate mediolateral, anteroposterior, and varus/valgus ligament function, and require augmentation and/or stem extensions due to inadequate bone stock;
- For use when both cruciate ligaments are excised and when greater varus/valgus constraint is required;
- ACCK can be used with augments and stem extensions, provides for personalized patient care by allowing a custom-fit of the femoral, tibial and patellar components independently, thus addressing the largest possible percentage of the population;
- High flexion up to 145 degree;
- Intraoperative flexibility with cylinder and offset extension stem;
- Ingenuity curvature design make sure any size of femoral implant fits with any size of tibial insert;
- Individual femoral and tibial augments with two different thicknesses;
- Bolt-Locking mechanism gives loosening zero chances.

Straight & Offset stems
with different diameter and
length options



Intraoperative flexibility
with two stem options:
straight and offset


Individual distal and
posterior femoral
anugments with two
thickness options

Tibial augment with
two thickness options



ACCK Femoral Condylar

Size Ref.	Cat. No.	Description	ML X AP (mm)
L1#	6915-1410	Name: ACCK Femoral Condylar	56 X 52
L2#	6915-1420		62 X 57
L3#	6915-1430		67 X 61
L4#	6915-1440		72 X 66
L5#	6915-1450		74 X 68
R1#	6915-1310	Regularly Supplied Material: Co-Cr-Mo	56 X 52
R2#	6915-1320		62 X 57
R3#	6915-1330		67 X 61
R4#	6915-1340		72 X 66
R5#	6915-1350		74 X 68



ACCK Femoral Condylar
(Co-Cr-Mo)


ACCK Tibial Tray

Size Ref.	Cat. No.	Description	ML X AP (mm)
A	7106-0002	Name: ACCK Tibial Tray Regularly Supplied Material: Co-Cr-Mo	59 X 38
B	7106-0003		63 X 41
C	7106-0004		67 X 43
D	7106-0005		71 X 46
E	7106-0006		75 X 48



ACCK Tibial Tray
(Co-Cr-Mo)

ACCK Tibial Insert



ACCK Tibial Insert
(UHMWPE)

Size Ref.	Cat. No.	Description	ML X AP (mm)	Tibial Tray	Width C (mm)	Height D (mm)
SA10mm	7224-1010	Name: ACCK Tibial Insert Regularly Supplied Material: UHMWPE	54 x 35	A	18.2	31
SA12mm	7224-1012		54 x 35		18.2	33
SA14mm	7224-1014		54 x 35		18.2	35
SA16mm	7224-1016		54 x 35		18.2	37
SA18mm	7224-1018		54 x 35		18.2	39
SA20mm	7224-1020		54 x 35		18.2	41
SA22mm	7224-1022		54 x 35		18.2	43
BC10mm	7224-1510		63 x 40	B, C	18.2	31
BC12mm	7224-1512		63 x 40		18.2	33
BC14mm	7224-1514		63 x 40		18.2	35
BC16mm	7224-1516		63 x 40		18.2	37
BC18mm	7224-1518		63 x 40		18.2	39
BC20mm	7224-1520		63 x 40		18.2	41
BC22mm	7224-1522		63 x 40		18.2	43
DE10mm	7224-2010		71 x 45	D, E	18.2	31
DE12mm	7224-2012		71 x 45		18.2	33
DE14mm	7224-2014		71 x 45		18.2	35
DE16mm	7224-2016		71 x 45		18.2	37
DE18mm	7224-2018		71 x 45		18.2	39
DE20mm	7224-2020		71 x 45		18.2	41
DE22mm	7224-2022		71 x 45		18.2	43